


SP.500*Spectrophotometer*

*for Ultraviolet and Visible
absorption measurements*

The frequency with which the Unicam SP.500 is mentioned in the literature is convincing proof of the sustained quality of its performance and of the wide variety of applications for which it is particularly suitable. For the maintenance of standards of both raw materials and the final product, for academic and medical research, wherever there is a problem of analysis or control, the easiest, most efficient, most economical answer may well be the Unicam SP.500 Spectrophotometer.

UNICAM

*An illustrated leaflet describing
the instrument in detail
will gladly be sent on request.*

UNICAM INSTRUMENTS LTD • ARBURY WORKS • CAMBRIDGE

U 132B

A new form...


... of the famous 62FM:
***Oertling announce
Model 62FMA***

- *ONE DIAL ONLY* for manipulation of weights up to 1 gram.
- *EXCEPTIONALLY QUICK* in operation.
- *SENSITIVITY* 0.2 mg. per division.

With these features it replaces its famous predecessor as the ideal "basic" precision balance for every type of laboratory. For full details of this and other new models please fill in the coupon on the opposite page.

OERTLING

PIONEERS IN BALANCE DESIGN FOR OVER 100 YEARS



Model 62FMA £89

(A)

*L. Oertling Ltd., Cray Valley Road, St. Mary Cray,
Orpington, Kent. Telephone: Orpington 5771 (5 lines).*

*Please send me full details of the latest range of Oertling Aperiodic and
Release-o-matic balances.*

Name

Address

.....

QUICKFIT apparatus off the shelf . . .

For the benefit of users of

‘QUICKFIT & QUARTZ’

Laboratory apparatus in London and the Home Counties we now carry a stock of over **7000** pieces of **QUICKFIT** apparatus covering a range of **420** items.

- Free and immediate delivery in our own vans.
-

Scientific Supplies Co. Ltd.

Scientific House : Vine Hill : London, E.C.1

Telephone: TERminus 2468

Over 36 years experience of complete Laboratory Furnishing.



If you use laboratory glassware —

YOU NEED THIS CATALOGUE

Enlarged and improved plant now enables us to give immediate delivery of the world's best laboratory glassware. A copy of our new catalogue will gladly be sent at your request.

QUICKFIT & QUARTZ LTD

STANDARD INTERCHANGEABLE LABORATORY
GLASSWARE

Orders and Enquiries to

Dept. O.A., "Quickfit" Works, Heart of STONE, Staffs

Phone: Stone 481

STANTON *Leadership*

The design and construction of every modern Precision Balance manifests the influence that STANTONS have exerted during the past few years.

STANTON were first with:

- Synthetic Sapphire (Corundum) Plates
- Full weight-loading balances in Britain
- Balances with special applications, e.g. Gas Density Determination, Thermo Recording work, Remote Control work.

STANTON now
present 2 new models


Model A.D.3.

Capacity: 200g.
Sensitivity: 0.1mg per $\frac{1}{2}$ division

This aperiodic projection reading balance provides for external weight loading of fractions up to 1g. by means of a single dial. The new graticule reads direct up to 100mg. resulting in greater weighing speed and less fatigue.



Fully descriptive illustrated folder on request



Model B.A.9.

Capacity: 200g.
Sensitivity: 0.1mg per $\frac{1}{2}$ division

This balance represents another contribution by STANTON towards greater balance room efficiency. Only four dials give weight loading up to 200gr.—Graticule reads direct up to 100mg. Complete weighing can thus be made in 20–30 seconds.

Stanton Balances are stocked by leading laboratory furnishers.



STANTON

PRECISION BALANCES


Stanton Instruments Ltd., 119 Oxford St., London, W.1. Gerrard 7533-4.

WHATMAN FILTER PAPER

GENERAL ANALYSIS: There is a large and varied choice of Whatman grades, including folded filter papers, single and double acid washed grades, and hardened filter papers designed to give great wet strength, resistance to alkalis and acids and low ash content.

CHROMATOGRAPHY: Whatman filter papers to be used in this field are specially selected and packed. In addition to the many grades available, Whatman Cellulose Powders in "Ashless" and "B" Quality, standard and coarse grades, offer another chromatographic medium for use in columns.

SPECIALITIES: For the convenience of users Whatman special aids for quick analysis include Drop Reaction Paper, Arsenic Test Papers. Ashless and Ordinary Clippings, Ashless Tablets, Soxhlet Extraction Thimbles.



Manufacturers :

W. & R. BALSTON LTD., MAIDSTONE, KENT

Sole Mill Representatives

H. REEVE ANGEL & CO. LTD.
9, Bridewell Place, London, E.C.4

A Photomultiplier

IN THE NEW SPEKKER FLUORIMETER provides greatly increased sensitivity in fluorescence measurement while retaining the special Spekker features of null point measurement and full compensation for light source variations. The null-point indicator is a robust spot galvanometer.



Write for full information to

HILGER & WATTS LTD

HILGER DIVISION

98, St. Pancras Way, Camden Road, London, N.W.1, England

Member of the Scientific Instrument Manufacturers' Association and of SCIEX

We have made available certain Analytical Reagents with ACTUAL BATCH ANALYSIS confirmed by INDEPENDENT Analysts of the highest standing: particulars of one example are given.

★ YOU ARE INVITED TO COMPARE THE PURITY WITH THAT GUARANTEED BY ANY COMPETING MAKER.

JUDACTAN ANALYTICAL REAGENT

TARTARIC ACID A.R.

(CHOH-COOH)₂

Mol. Wt. 150.09

ACTUAL BATCH ANALYSIS

(Not merely maximum impurity values)

Batch No. 23001

Ash	0.02%
Chloride (Cl)	0.0001%
Sulphate (SO ₄)	0.001%
Lead (Pb)	0.00002%
Iron (Fe)	0.00004%
Arsenic (As ₂ O ₃)	no reaction

The above analysis is based on the results, not of our own Control Laboratories alone, but also on the confirmatory Analytical Certificate issued by independent Consultants of international repute.

The General Chemical & Pharmaceutical Co., Ltd.
Chemical Manufacturers, Judex Works, Sudbury, Middlesex